

## Product datasheet for **AM31884PU-N**

### CD11a / ITGAL Rat Monoclonal Antibody [Clone ID: IBL-6/2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	IBL-6/2
Applications:	FC, IHC, IP
Recommended Dilution:	<b>Immunohistochemistry on frozen sections.</b> <b>Immunoprecipitation.</b> <b>Flow Cytometry.</b>
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	EL4 (mouse thymoma cell line) <u>Donor:</u> Wistar spleen <u>Fusion Partner:</u> Sp-2/0Ag14
Specificity:	This antibody reacts to CD11a / ITGAL.
Formulation:	PBS State: Purified State: Liquid purified Ig fraction Stabilizer: EIA grade BSA Preservative: 0.02% sodium azide
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">P24063</a>
Background:	The CD11a is expressed on leukocytes and it mediates cell-cell and cell-adhesion by binding to CD54 and ICAM-2.
Synonyms:	Integrin alpha-L, LFA1, LFA-1



[View online »](#)