

## Product datasheet for **AM31885BT-N**

### Havcr2 Rat Monoclonal Antibody [Clone ID: 25F.1D6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	25F.1D6
Applications:	FC
Recommended Dilution:	Immunohistochemistry. Flow Cytometry.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	Anti-mouse Tim-3 (T cell immunoglobulin domain, mucin domain 3) monoclonal antibody reacts with mouse Tim-3.
Formulation:	PBS containing 0.02% sodium azide (NaN <sub>3</sub> ) as preservative and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/mL. Label: Biotin State: Liquid purified IgG fraction
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	hepatitis A virus cellular receptor 2
Database Link:	<a href="#">Entrez Gene 171285 Mouse</a> <a href="#">Q8VIM0</a>



[View online »](#)

**Background:** This 281 amino acid type 1 transmembrane protein contains an immunoglobulin and mucin-like domain in its extra-cellular portion and a tyrosine phosphorylation motif in its cytoplasmic portion. Tim-3 is a member of the immunoglobulin super-family and is expressed on Th1 lymphocytes and CD11b+ macrophages. Tim-3 is a T helper type 1 (TH1)-specific molecule that regulates TH1 mediated auto and alloimmune responses and the induction of peripheral immunological tolerance.

**Synonyms:** TIMD-3, TIM3, TIM-3, T-cell membrane protein 3, HAVCR-2