

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for AM31885FC-N

Havcr2 Rat Monoclonal Antibody [Clone ID: 25F.1D6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	25F.1D6
Applications:	FC, IHC
Recommended Dilution:	Immunohistochemistry. Flow Cytometry.
Reactivity:	Mouse
Host:	Rat
lsotype:	lgG2a
Clonality:	Monoclonal
Specificity:	Anti-mouse Tim-3 (T cell immunoglobulin domain, mucin domain 3) monoclonal antibody reacts with mouse Tim-3.
Formulation:	PBS containing 0.02% sodium azide (NaN39 as preservative and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/mL. Label: FITC State: Liquid purified IgG fraction Label: Fluorescein isothiocyanate isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	hepatitis A virus cellular receptor 2
Database Link:	<u>Entrez Gene 171285 Mouse</u> <u>Q8VIM0</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Havcr2 Rat Monoclonal Antibody [Clone ID: 25F.1D6] – AM31885FC-N
Background:	This 281 amino acid type 1 transmembrane protein contains an immunoglobulin and mucin- like domain in its extra-cellular portion and a tyrosine phosphorylation motif in its cytoplasmic portion. Tim-3 is a member of the immunoglobulin super-family and is expressed on Th1 lymphocytes and CD11b+ macrophages. Tim-3 is a T helper type 1 (TH1)-specific molecule that regulates TH1 mediated auto and alloimmune responses and the induction of peripheral immunological tolerance.
Synonyms:	TIMD-3, TIM3, TIM-3, T-cell membrane protein 3, HAVCR-2

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US