

Product datasheet for AM31928PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Neuropeptide Y (NPY) (29-97) Mouse Monoclonal Antibody [Clone ID: 2C10]

Product data:

Product Type: Primary Antibodies

Clone Name: 2C10

Applications: ELISA, IHC, WB

Recommended Dilution: ELISA.

Immunohistochemistry on Paraffin Sections: 5 µg/ml.

Western Blot: 1/500 - 1/1000.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant protein encoding amino acids 29-97 of human NPY

Specificity: This antibody recognizes Human NPY / Neuropeptide Y at aa 29-97.

Formulation: PBS, pH 7.2

State: Purified

State: Liquid purified Ig fraction

Concentration: lot specific

Purification: Protein A Chromatography

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: neuropeptide Y

Database Link: Entrez Gene 4852 Human

P01303

Background: Neuropeptide Y (NPY) is a peptide belonging to the PP-family and occurs in neurons and

adrenal medullary cells. The brain contains large quantities of NPY and it is also found in the peripheral nervous system, where it coexists with noradrenaline in sympathetic fibers. NPY

inhibits gut motility and causes vasoconstriction. Pheochromocytomas contain NPY.






Synonyms: OTTHUMP00000201946; OTTHUMP00000201947; PYY4


Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Adipocytokine signaling pathway

Product images:




Western Blot analysis of NPY expression in transfected 293T cell line by NPY monoclonal antibody, clone 2C10.



Detection limit for recombinant GST tagged NPY is approximately 3ng/ml as a capture antibody.





Human Adrenal Medulla: Formalin-Fixed, Paraffin-Embedded (FFPE)