

Product datasheet for AM31967PU-N

Il10 Mouse Monoclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: In a Sandwich ELISA (assuming 100 µl/well), a concentration of 8.0-16.0 µg/ml of this antibody will detect at least 500 pg/well of recombinant Rat IL-10 when used in conjunction with Biotinylated antigen affinity purified anti-Rat IL-10 (<i>CatNo</i> PP1115B1/PP1115B2) as the detection antibody at a concentration of at least 0.5-1.0 µg/ml. Western Blot: To detect Rat IL-10 by Western Blot analysis this antibody can be used at a concentration of 0.5-1.0 µg/ml. When used in conjunction with compatible secondary reagents the detection limit for recombinant Rat IL-10 is 0.5-1.0 µg/lane, under reducing or non-reducing conditions.
Reactivity:	Rat
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	<i>E.coli</i> derived Recombinant Rat IL-10 (<i>CatNo</i> PA041).
Specificity:	This antibody recognizes Interleukin-10 / IL10
Formulation:	PBS State: Aff - Purified State: Lyophilized purified tissue culture
Reconstitution Method:	Centrifuge the vial prior to opening and restore in water to a concentration of 1.0 mg/ml.
Purification:	Affinity Chromatography Protein A
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 10
Database Link:	<u>Entrez Gene 25325 Rat</u> <u>P29456</u>




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Background:IL-10 is an immunosuppressive cytokine produced by a variety of mammalian cell types
including macrophages, monocytes, T cells, B cells and keratinocytes. IL-10 inhibits the
expression of proinflammatory cytokines such as IL-1 and TNF-alpha. Like IL-4, IL-10
enhances humoral immune responses and attenuates cell-mediated immune reactions.
Human IL-10 is active on murine cells, but murine IL-10 is inactive on human cells.

Synonyms: IL-10, CSIF, TGIF

Product images:



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US