

## Product datasheet for **AM31975PU-N**


### Rhotekin 2 (RTKN2) (1-164) Mouse Monoclonal Antibody [Clone ID: 2C2]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	2C2
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry on Paraffin Sections:</b> 1/50. <b>Western Blot.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	RTKN2 (AAH25765, 1 a.a. ~164 a.a) full-length recombinant protein with GST tag
Specificity:	Recognizes RTKN2 at aa 1-164.
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	rhotekin 2
Database Link:	<a href="#">Entrez Gene 219790 Human</a> <a href="#">Q8IZC4</a>
Synonyms:	PLEKHK1



[View online »](#)

**Product images:**

Detection limit for recombinant GST tagged RTKN2 is 0.3 ng/ml as a capture antibody.



Human Pancreas: Formalin-Fixed, Paraffin-Embedded (FFPE)