

Product datasheet for **AM31997PU-N**

Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: MF4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: MF4]
Clone Name:	MF4
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Chicken
Host:	Mouse
Immunogen:	Polyclonal IgY (Promega, REF G1161).
Isotype:	IgG1
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgY fraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or Add 0.05% Sodium Azide for longer storage.



[View online »](#)