

Product datasheet for **AM31998PU-N**

Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 08C]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 08C]
Clone Name:	08C
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Immunoblotting.
Reactivity:	Chicken
Host:	Mouse
Immunogen:	Polyclonal IgY (Promega, REF G1161).
Isotype:	IgG1
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgY fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or Add 0.05% Sodium Azide for longer storage.



[View online »](#)

Product images:




Figure 1: Specificity of Anti- chicken IgG (IgY) Heavy Chain Immunoglobulin, clone 0.8C, determined by Immunoblot analysis. Blot contains Na-alginate precipitated yolk IgY (predicted band size [heavy chain IgY] 70 kDa). Incubated with antibody fraction (culture supernatant) 10X diluted in PBS containing 0.05% tween-20 and 5% non fat dry milk.