

## Product datasheet for **AM32000PU-N**

### Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 9A6]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 9A6]
Clone Name:	9A6
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunoblotting.</b>
Reactivity:	Chicken
Host:	Mouse
Immunogen:	Polyclonal IgY (Promega, REF G1161).
Isotype:	IgG2a
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgY fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or Add 0.05% Sodium Azide for longer storage.



[View online »](#)

## Product images:




Figure 1: Specificity of Anti- chicken IgG (IgY) Heavy Chain Immunoglobulin, clone 9.A6, determined by immunoblot analysis. Blot contains Na-alginate precipitated yolk IgY (predicted band size [heavy chain IgY] 70 kDa). Incubated with antibody fraction (culture supernatant) 10X diluted in PBS containing 0.05% tween-20 and 5% non fat dry milk.