

Product datasheet for **AM32009PU-N**

Collagenase G / ColG Mouse Monoclonal Antibody [Clone ID: 1D2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1D2
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/250-1/500 (guidelines dilution). Western Blot: 1/2000-1/5000.
Reactivity:	Clostridium histolyticum
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Collagenase from <i>Clostridium histolyticum</i>
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: One year from despatch.
Database Link:	Q9X721
Background:	This Collagenase is typically used with a neutral protease for tissue disruption.



[View online »](#)