

Product datasheet for **AM32132BT-N**

IP10 (CXCL10) Mouse Monoclonal Antibody [Clone ID: 6D4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	6D4
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	Label: Biotin State: Purified Storage Buffer: PBS with 0.02% Sodium Azide and 0.1% BSA
Conjugation:	Biotin
Gene Name:	C-X-C motif chemokine ligand 10
Database Link:	Entrez Gene 3627 Human P02778
Synonyms:	C-X-C motif chemokine 10, Small-inducible cytokine B10, Gamma-IP10, INP10, SCYB10



[View online »](#)