

## Product datasheet for **AM32194FC-N**

### Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]
Clone Name:	LO-MG-7
Applications:	ELISA, IF, IHC
Recommended Dilution:	<b>Immunoassay.</b> <b>ELISA:</b> good binding on plastic. <b>Immunohistochemistry on frozen sections.</b> <b>Indirect Immunofluorescence.</b>
Reactivity:	Mouse
Host:	Rat
Isotype:	IgM
Formulation:	PBS Label: FITC State: Purified State: Liquid purified Ig fraction Preservative: 0.09% sodium azide
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)