

Product datasheet for **AM32194HR-N**

Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]
Clone Name:	LO-MG-7
Applications:	ELISA, IHC
Recommended Dilution:	Immunoassay. ELISA: good binding on plastic. Immunohistochemistry on Frozen Sections.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgM
Formulation:	PBS Label: HRP State: Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol
Concentration:	lot specific
Conjugation:	HRP
Storage:	Store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)