

Product datasheet for **AM32194PU-N**

Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG-7]
Clone Name:	LO-MG-7
Applications:	ELISA, IHC
Recommended Dilution:	ELISA (good binding on plastic) LO-MG-7 purified can be used at 0.6 µg/ml to obtain an OD of more than 1 on Mouse IgG2a en IgG2b. Immunoassay: Useful as Secondary antibody. Immunohistochemistry on Frozen Sections.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgM
Formulation:	PBS State: Purified State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)