

Product datasheet for **AM32195BT-N**

Mouse IgA (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MA-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MA-7]
Clone Name:	LO-MA-7
Applications:	ELISA, IHC
Recommended Dilution:	Immunoassay. ELISA: good binding on plastic. Immunohistochemistry on frozen sections.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgM
Formulation:	PBS Label: Biotin State: Purified State: Liquid purified Ig fraction Stabilizer: 50% glycerol Preservative: 0,09% sodium azide
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)