

Product datasheet for **AM32197PU-N**

Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 7C2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 7C2]
Clone Name:	7C2
Applications:	ELISA, IP, WB
Recommended Dilution:	Immunoassay. Western Blot (non-reduced). The typical starting working dilution is 1/10. Immunoprecipitation.
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified Ig fraction (0.2 µm filtered) Stabilizer: 0,1% BSA Preservative: 0,02% sodium azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)