

Product datasheet for **AM32203PU-N**

Biotin Mouse Monoclonal Antibody [Clone ID: BIO-8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BIO-8
Applications:	IHC
Recommended Dilution:	Immunohistochemistry.
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Specificity:	The antibody recognizes an epitope that is exposed by protein-biotin and nucleic acid-biotin conjugates.
Formulation:	State: Purified State: Liquid tissue culture supernatant Preservative: 0,09% sodium azide
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Vitamin B7, Vitamin H



[View online »](#)