

Product datasheet for AM32216FC-N

Human IgM Mouse Monoclonal Antibody [Clone ID: ICL-931]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Mouse Monoclonal Antibody [Clone ID: ICL-931]
Clone Name:	ICL-931
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on Frozen Sections: 1/40. Staining of Human cells and tissue sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgM.
lsotype:	lgG1
Formulation:	0.01M Sodium Phosphate, pH 7.5 with 50% Glycerol/0.1M NaCl Label: FITC State: Liquid purified Ig fraction Label: Fluorescein Isothiocyanate Isomer I
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C.

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US