

Product datasheet for **AM32216PU-N**

Human IgM Mouse Monoclonal Antibody [Clone ID: ICL-931]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Mouse Monoclonal Antibody [Clone ID: ICL-931]
Clone Name:	ICL-931
Applications:	ELISA, IHC, R
Recommended Dilution:	Dot Blot. ELISA. Radioimmunoassay. Immunohistochemistry on frozen sections: 1/40.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgM
Isotype:	IgG1
Formulation:	1% mannitol 0.15 M NaCl, 0.01 M Naphosphate at pH 7.5 State: Purified State: Lyophilized purified Ascites
Reconstitution Method:	Reconstitute with distilled water.
Purification:	Sulphate precipitation and DEAE-cellulose chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)