

Product datasheet for **AM32217FC-N**

Human IgG Mouse Monoclonal Antibody [Clone ID: EFE-565]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Mouse Monoclonal Antibody [Clone ID: EFE-565]
Clone Name:	EFE-565
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/40.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgG
Isotype:	IgG1
Formulation:	50% glycerol/0.1 M NaCl, 0.01 M Na-phosphate at pH 7.5 Label: FITC State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)