

Product datasheet for **AM32217HR-N**

Human IgG Mouse Monoclonal Antibody [Clone ID: EFE-565]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Mouse Monoclonal Antibody [Clone ID: EFE-565]
Clone Name:	EFE-565
Applications:	ELISA, IHC, R
Recommended Dilution:	Dot Blot. ELISA. Radioimmunoassay. Immunohistochemistry on Frozen Sections: 1/40.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified Human IgG
Isotype:	IgG1
Formulation:	50% Glycerol 0.15M NaCl, 0.01M Sodium Phosphate at pH 7.5 Label: HRP State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Gel-filtration and ion-exchange chromatography to remove peroxidase and unlabeled antibodies
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C.



[View online »](#)