

## Product datasheet for **AM32224FC-N**

### Human IgA Mouse Monoclonal Antibody [Clone ID: BFG-267]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA Mouse Monoclonal Antibody [Clone ID: BFG-267]
Clone Name:	BFG-267
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/40.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified $\alpha$ -chain (Fc) of human IgA.
Isotype:	IgG1
Formulation:	50% glycerol/0.1 M NaCl, 0.01 M Na-phosphate at pH 7.5 Label: FITC State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)