

Product datasheet for **AM32224PU-N**

Human IgA Mouse Monoclonal Antibody [Clone ID: BFG-267]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA Mouse Monoclonal Antibody [Clone ID: BFG-267]
Clone Name:	BFG-267
Applications:	ELISA, IHC, R
Recommended Dilution:	ELISA. Radioimmunoassay. Dot Blot. Immunohistochemistry on frozen sections: 1/40.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified α -chain (Fc) of human IgA.
Isotype:	IgG1
Formulation:	1% mannitol 0.15 M NaCl, 0.01 M Naphosphate at pH 7.2 State: Purified State: Lyophilized purified Ascites
Reconstitution Method:	Reconstitute with distilled water.
Purification:	Sulphate precipitation and DEAE-cellulose chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody at 2-8°C.



[View online »](#)