

Product datasheet for **AM32284PU-N**

ACTH (POMC) Mouse Monoclonal Antibody [Clone ID: 2F6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	2F6
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on frozen and paraffin sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	A BALB/c mouse was immunized with ACTH synthetic peptide (1-24 fragment) conjugated KLH. Splenocytes were fused with SP-2/0.
Specificity:	The mAb 2F6 demonstrates broad species cross reactivity. It localizes ACTH in immunohistochemical techniques on formalin-fixed, paraffin-embedded or frozen sections
Formulation:	PBS State: Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	proopiomelanocortin
Database Link:	Entrez Gene 5443 Human P01189
Synonyms:	Pro-opiomelanocortin, Corticotropin-Lipotropin



[View online »](#)