

Product datasheet for **AM32407PU-N**

HLAB (HLA-B) Mouse Monoclonal Antibody [Clone ID: JOAN-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	JOAN-1
Applications:	IHC, WB
Recommended Dilution:	Western Blot. Immunohistochemistry on frozen sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	NS1 mouse myeloma cells were fused with splenocytes from a BALB/c mouse immunized intraperitoneally with normal human peripheral blood lymphocytes.
Specificity:	This antibody recognizes an intralocus determinant present on a limited number of HLA-B locus encoded gene products.
Formulation:	PBS State: Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	major histocompatibility complex, class I, B
Database Link:	Entrez Gene 3106 Human P01889
Synonyms:	HLAB, MHC class I antigen B7



[View online »](#)