

## Product datasheet for **AM32410SU-N**

### CD3E Mouse Monoclonal Antibody [Clone ID: SPV-T3b]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	SPV-T3b
Applications:	IF, IHC
Recommended Dilution:	<b>Immunofluorescence.</b> <b>Immunohistochemistry on frozen sections.</b> Recommended dilution: 1/10 - 1/20.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	HG-38
Specificity:	This antibody reacts with human CD3, 19-29 kD, present on all resting and activated human T-cells, on Tleukemia cells and a proportion of human thymocytes.
Formulation:	PBS State: Ascites State: Liquid purified Ascites Preservative: 0,09% sodium azide
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD3e molecule
Database Link:	<a href="#">Entrez Gene 916 Human</a> <a href="#">P07766</a>
Synonyms:	T3/Leu-4



[View online »](#)