

Product datasheet for **AM32411SU-N**

CD8A Mouse Monoclonal Antibody [Clone ID: FK18]

Product data:

Product Type:	Primary Antibodies
Clone Name:	FK18
Applications:	IF
Recommended Dilution:	Immunofluorescence: 1/10. Preferably in PBS plus 1%BSA and 0.02% sodium azide is advised
Reactivity:	Human, Monkey
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Specificity:	This antibody reacts with CD8, the 33 kD T8 antigen present on the human cytotoxic/suppressor T-cell subset. It strongly blocks the cytolytic activity of T8-positive CTLs. It also reacts with a subpopulation of Rhesus monkey lymphocytes which is identical to the population of human T8 positive cells.
Formulation:	liquid antibody solution with sodium azide.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD8a molecule
Database Link:	Entrez Gene 925 Human P01732
Synonyms:	CD8 alpha chain, CD8A, MAL



[View online »](#)