

Product datasheet for **AM32474PU-N**

Plasmodium falciparum Histidine Rich Protein-2 (HRP2) Mouse Monoclonal Antibody [Clone ID: B1242M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B1242M
Applications:	ELISA, LF
Recommended Dilution:	ELISA. Lateral Flow. <i>Recommended antibody pairs for Sandwich Immunoassay.</i> Use AM32474PU-N (Clone B1242M) as Capture Antibody and AM32475PU-N (Clone B1245M) or AM32476PU-N (Clone B1243M) as Detection Antibody.
Reactivity:	Plasmodium falciparum
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native <i>Plasmodium falciparum</i> HRP-2.
Specificity:	This Monoclonal antibody (B1242M) reacts with HRP-2 from <i>P. falciparum</i> and Does not react with other <i>Plasmodium</i> species.
Formulation:	0.01M PBS, pH 7.3 State: Purified State: Liquid purified IgG fraction Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.



[View online »](#)