

## Product datasheet for **AM32475PU-N**

### Plasmodium falciparum Histidine Rich Protein-2 (HRP2) Mouse Monoclonal Antibody [Clone ID: B1245M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B1245M
Applications:	ELISA, LF
Recommended Dilution:	<b>ELISA.</b> <b>Lateral Flow.</b> <i>Recommended antibody pairs for Sandwich Immunoassay.</i> Use AM32474PU-N (Clone B1242M) as Capture Antibody and AM32475PU-N (Clone B1245M) as Detection Antibody.
Reactivity:	Plasmodium falciparum
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native <i>Plasmodium falciparum</i> HRP-2
Specificity:	This Monoclonal antibody (Clone B1245M) reacts with HRP-2 from <i>P. falciparum</i> and <b>Does not</b> react with other <i>Plasmodium</i> species.
Formulation:	0.01M PBS, pH 7.3 State: Purified State: Liquid purified IgG fraction Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.



[View online »](#)