

## Product datasheet for AM33351PU-S

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## TNF alpha (TNF) Mouse Monoclonal Antibody [Clone ID: 4C6-H8]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 4C6-H8

**Applications:** FC, IF, IHC, IP

**Recommended Dilution: ELISA:** Use free BSA Antibody for coaating.

**Flow Cytometry:** 0.5-1 µg/million cells. **Immunofluorescence:** 0.5-1 µg/ml.

**Immunoprecipitation (Not tested):** 1-2 μg/500 μg protein lysate.

Immunohistochemistry on Frozen and Formalin-Fixed Sections: 2-4 µg/ml for 30 minutes

at RT.

Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris Buffer with

1mM EDTA, pH 9.0, for 10-20 min followed by cooling at RT for 20 minutes.

**Positive Control**: HeLa, HL-60, or A431 cells. Macrophages in lymph node or Tonsil.

**Reactivity:** Canine, Feline, Human, Mouse, Rabbit, Rat, Zebrafish

**Host:** Mouse

Isotype: IgM

Clonality: Monoclonal

**Immunogen:** A hexadecapeptide corresponding to amino acids 115-130 of Human TNF-alpha, conjugated

to Thyroglobulin.

**Specificity:** Recognizes TNF-alpha.

Clone 4C6-H8 reacts with Tumor Necrosis Factor-alpha. It reacts on paraffin sections with macrophages, in which the cytoplasma is stained. Some keratinocytes are also positive

(tonsils).

**Cellular Localization**: Cytoplasmic and extracellular (secreted).

Formulation: 10mM PBS

State: Purified

State: Liquid purified IgG fraction from Bioreactor Concentrate

Stabilizer: 0.05% BSA

Preservative: 0.05% Sodium Azide

**Concentration:** lot specific

**Purification:** Protein A/G Chromatography





## TNF alpha (TNF) Mouse Monoclonal Antibody [Clone ID: 4C6-H8] - AM33351PU-S

Conjugation: Unconjugated

**Storage:** Store undiluted at 2-8°C.

**Stability:** Shelf life: one year from despatch.

**Predicted Protein Size:** 17 kDa

Gene Name: tumor necrosis factor

Database Link: Entrez Gene 7124 Human

P01375

**Background:** Tumor Necrosis Factor Alpha (TNF alpha) is a protein secreted by lipopolysaccharide-

stimulated macrophages, and causes tumor necrosis when injected into tumor bearing mice.

TNF alpha is believed to mediate pathogenic shock and tissue injury associated with endotoxemia. TNF alpha exists as a multimer of two, three, or five non-covalently linked units, but shows a single 17kDa band following SDS PAGE under non-reducing conditions. TNF

alpha is closely related to the 25kDa protein Tumor Necrosis Factor beta (lymphotoxin), sharing the same receptors and cellular actions. TNF alpha causes cytolysis of certain

transformed cells, being synergistic with interferon gamma in its cytotoxicity. Although it has little effect on many cultured normal human cells, TNF alpha appears to be directly toxic to vascular endothelial cells. Other actions of TNF alpha include stimulating growth of human fibroblasts and other cell lines, activating polymorphonuclear neutrophils and osteoclasts, and induction of interleukin 1, prostaglandin E2 and collagenase production. TNF alpha is currently being evaluated in treatment of certain cancers and AIDS Related Complex.

**Synonyms:** TNF, TNF-a, TNFA, TNFSF2, Cachectin