

## Product datasheet for **AM33379AF-N**

### galectin 9 (LGALS9) Mouse Monoclonal Antibody [Clone ID: B-K38]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-K38
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant Human Galectin-9. <b>Hybridoma:</b> Myeloma X63/AG.8653 x Balb/c node cells.
Specificity:	This Monoclonal Antibody recognizes both natural and recombinant Human Galectin-9 protein. Other species not tested.
Formulation:	PBS, Sterile-filtered through 0.22 $\mu$ m. Carrier and preservative free. State: Azide Free State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Stability:	Shelf life: one year from despatch.
Gene Name:	galectin 9
Database Link:	<a href="#">Entrez Gene 3965 Human</a> <a href="#">O00182</a>
Synonyms:	LGALS9, Galectin9, HOM-HD-21, Ecalectin



[View online »](#)