

## Product datasheet for **AM33403PU-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: B1521M]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: B1521M]
Clone Name:	B1521M
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	IgG from Human Serum.
Isotype:	IgG
Formulation:	PBS, pH 7.4 State: Purified State: Liquid purified IgG fraction from Cell Culture (> 90% pure by SDS-PAGE) Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)