

Product datasheet for **AM33404PU-N**

Human IgM (Mu chain, Fab-Region) Mouse Monoclonal Antibody [Clone ID: B1545M]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain, Fab-Region) Mouse Monoclonal Antibody [Clone ID: B1545M]
Clone Name:	B1545M
Applications:	ELISA
Recommended Dilution:	This antibody can pair with itself in ELISA, EIA . Can be used for specific IgM antibodies detection and for total IgM level detection. Can also be used as a capture antibody for total IgM level detection.
Reactivity:	Human
Host:	Mouse
Immunogen:	IgM from Human Serum.
Isotype:	IgG2b
Formulation:	PBS, pH 7.4 State: Purified State: Liquid purified IgG fraction from Cell Culture (> 90% pure by SDS-PAGE) Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)