

Product datasheet for **AM33468PU-N**

Herpes simplex Virus 1 / HSV1 Mouse Monoclonal Antibody [Clone ID: B323M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B323M
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human Alphaherpesvirus 1
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	HSV1 (strain HFEM).
Specificity:	The antibody recognizes HSV1.
Formulation:	PBS, pH 7.5 State: Purified State: Liquid purified Ig fraction from Culture (> 90% pure by SDS-PAGE) Preservative: 15mM Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Synonyms:	HSV-1, HHV1, HHV-1, Human Herpes Virus 1



[View online »](#)