

## Product datasheet for **AM33470PU-N**

### Herpes simplex Virus 2 / HSV2 Mouse Monoclonal Antibody [Clone ID: B318M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B318M
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human Alphaherpesvirus 2
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	HSVII (strain 25766).
Specificity:	The antibody recognizes HSV2.
Formulation:	PBS, pH 7.5 State: Purified State: Liquid purified Ig fraction from Culture (> 90% pure by SDS-PAGE) Preservative: 15mM Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Stability:	Shelf life: one year from despatch.
Synonyms:	HSV-2, HHV2, HHV-2, Human Herpes Virus 2



[View online »](#)