

Product datasheet for **AM33472PU-N**

Influenza B (Nucleoprotein) Mouse Monoclonal Antibody [Clone ID: B312M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B312M
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Influenza B Virus
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Influenza B/Lee/40 and Influenza B/Singapore/-222/79 virus.
Specificity:	This antibody recognizes Influenza B Virus Nucleoprotein (species specific conserved epitope).
Formulation:	PBS, pH 7.5 State: Purified State: Liquid purified Ig fraction (> 90% pure by SDS-PAGE) Preservative: 15 mM Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza B Virus



[View online »](#)