

Product datasheet for **AM33479AF-N**

IL21 Mouse Monoclonal Antibody [Clone ID: B-P41]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-P41
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant Human IL-21 Hybridoma: Myeloma X63/AG.8653 x Balb/c spleen cells.
Specificity:	This Monoclonal Antibody recognizes both recombinant and natural Human IL-21. Other species not tested.
Formulation:	PBS, Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 21
Database Link:	Entrez Gene 59067 Human Q9HBE4
Synonyms:	IL-21, Za11



[View online »](#)