

## Product datasheet for **AM33481AF-N**

### IL32 Mouse Monoclonal Antibody [Clone ID: B-B52]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-B52
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant Human IL-32 gamma. <b>Hybridoma:</b> Myeloma X63/AG.8653 x Balb/c spleen cells.
Specificity:	This Monoclonal Antibody recognizes both recombinant and natural Human IL-32 32 (all isoforms). Other species not tested.
Formulation:	PBS, Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 32
Database Link:	<a href="#">Entrez Gene 9235 Human P24001</a>
Synonyms:	IL-32, NK4, TAIF



[View online »](#)