

## Product datasheet for **AM39045SU-N**

### LOC105369261 Mouse Monoclonal Antibody [Clone ID: HIS17]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS17
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	The antibody reacts with an epitope of CD43 (Leukosialin) on the surface of plasma cells, granulocytes, monocytes, macrophages, natural killer cells, thymocytes and most peripheral T cells. CD43 is not expressed on resting peripheral B cells.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Database Link:	<a href="#">P16150</a>
Synonyms:	Leukocyte sialoglycoprotein, Sialophorin, Galactoglycoprotein, SPN



[View online »](#)