

## Product datasheet for **AM39048SU-N**

### MHC Class II (non polymorphic) Mouse Monoclonal Antibody [Clone ID: HIS28]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS28
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	The antibody recognizes non-polymorphic MHC Class II antigen, all haplotypes.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	MHC Class 2, HLA Class II



[View online »](#)