

Product datasheet for **AM39049SU-N**

MHC Class II RT1B Mouse Monoclonal Antibody [Clone ID: HIS19]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS19
Applications:	IF, IHC, IP
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections. Immunohistochemistry on Paraffin and Resin sections: This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase. Immunoprecipitation.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Peyer's patch lymphocytes from AO strain rats.
Specificity:	The antibody recognizes polymorphic MHC Class II (RT1B), all haplotypes except RT1n.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	HLA Class II, MHC Class 2 RT1B



[View online »](#)