

## Product datasheet for **AM39050SU-N**

### **B cells and T cells (Pan) Mouse Monoclonal Antibody [Clone ID: HIS23]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	HIS23
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	The antibody recognizes rat pan B cells and pan T cells. Other species not tested.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)