

Product datasheet for **AM39052SU-N**

LAMA tumor / HIS30 Mouse Monoclonal Antibody [Clone ID: HIS30]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS30
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	The antibody reacts with LAMA pre-B cell leukemia, RStC tumor and Roser leukemia cells. Reactivity with normal cells: a subset of adrenal cortex (Zona glomerulosa).
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	LAMA pre-B cell tumor



[View online »](#)