

## Product datasheet for **AM39053SU-N**

### LAMA tumor subset / HIS31 Mouse Monoclonal Antibody [Clone ID: HIS31]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS31
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	The antibody recognizes a LAMA pre-B cell tumor subset.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Sore undiluted at 2-8°C.
Synonyms:	LAMA pre-B cell tumor



[View online »](#)