

Product datasheet for **AM39056SU-N**

MHC Class II (like thymus) Mouse Monoclonal Antibody [Clone ID: HIS34]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS34
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	The antibody recognizes MHC Class II antigen - like thymus, not peripheral B cells.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	MHC Class 2, HLA Class II



[View online »](#)