

Product datasheet for **AM39058SU-N**

Macrophages (ED2-like antigen) Mouse Monoclonal Antibody [Clone ID: HIS36]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS36
Applications:	IF, IHC
Recommended Dilution:	Flow Cytometry: not tested Immunofluorescence. Immunohistochemistry on Frozen tissue Sctions (1/20).
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	The antibody recognizes rat macrophages by reacting with an ED2-like antigen. Other species not tested.
Formulation:	State: Supernatant State: Liquid tissue culture Supernatant containing 8% (v/v) FBS Preservative: 0.09% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2–8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Macrophage marker

