

Product datasheet for **AM39059SU-N**

Thymic cortex epithelium Mouse Monoclonal Antibody [Clone ID: HIS37]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS37
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	The antibody recognizes rat thymic cortex epithelium. Other species not tested.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	Thymic epithelial cells



[View online »](#)