

Product datasheet for **AM39061SU-N**

Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: HIS40]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: HIS40]
Clone Name:	HIS40
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen tissue sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)