

Product datasheet for **AM39066SU-N**

T cells / Thymocytes (QCA-1) Mouse Monoclonal Antibody [Clone ID: HIS45]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HIS45
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Subset T cells (resting T cells), subset thymocytes. QCA-1 antigen.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)