

## Product datasheet for **AM39074SU-N**

### **B-Cells (corona) Mouse Monoclonal Antibody [Clone ID: HIS59]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	HIS59
Applications:	IF, IHC
Recommended Dilution:	Flow cytometry. Immunofluorescence. Immunohistochemistry on frozen sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	The antibody recognizes rat corona B-cells. Other species not tested.
Formulation:	State: Supernatant State: Liquid cell culture supernatant containing 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Synonyms:	B-cell marker, B cells



[View online »](#)